



Imgen BioSciences, Inc.

151 Martine Street, Fall River, MA 02723, USA Tel: 1 (781) 850-2047 www.imgenbiosciences.com

A-2101-2: Rat Anti-mouse IgG1 Monoclonal Antibody Clone KT91

Product Name:	Rat Anti- mouse IgG1 Monoclonal Antibody Clone KT91
Catalogue No:	A-2101-2
Immunogen:	Mouse IgG1
Source/Host:	Rat
Purity/Purification:	Protein G purified
Clone:	Monoclonal, Clone KT91
Antibody Class:	IgG
Species reactivity:	Not cross-reactive with mouse IgG2a, IgG2b, IgG3 and IgM
Form:	Liquid
Concentration:	1.0 mg/ml; 0.1% Sodium Azide
Size:	0.2 mg
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Applications:	ELISA
Shipping:	May be shipped with ice packs or dry ice.

Brief description about the rat anti-mouse IgG1 monoclonal antibody clone KT91:

This monoclonal antibody was shown not cross-reactive with mouse IgG2a, IgG2b, IgG3 and IgM by ELISA.

(Produced by Imgen BioSciences, Inc., April, 2013)