



Imgen BioSciences, Inc.

151 Martine Street, Fall River, MA 02723, USA Tel: 1 (781) 850-2047 www.imgenbiosciences.com

A-2111-1: Goat Anti-Mouse IgG1 Antibody

Product Name:	Goat Anti-Mouse IgG1 Antibody
Catalogue No:	A-2111-1
Immunogen:	Highly purified mouse IgG1
Source/Host:	Goat
Purity/Purification:	Sequential absorption with IgG2a, IgG2b and IgG3 immunoabsorbants until it is IgG1 specific followed by immunoaffinity purification with mouse IgG1.
Clone:	Polyclonal
Species reactivity:	IgG1
Form:	Liquid
Concentration:	1.0 mg/ml (in PBS containing 0.1% sodium azide)
Size:	0.2 mg
Storage:	Store at 2-8 °C
Applications:	ELISA, WB: Tested
Shipping:	May be shipped with ice packs or dry ice.

(Produced by Imgen BioSciences, Inc. April 2013)